



dutch cases

# DISCLAIMER

---

Please be aware that different healthcare providers may have varying approaches to lab testing and interpretation. The selection of specific tests, methodologies, and treatment recommendations can differ based on the provider's training, experience, and the individual needs of the patient. Providers should always use their best clinical judgment when making decisions for patient care.

The medical information in this case study is provided as an information resource only and is not to be used or relied on for any diagnostic or treatment purposes. This lecture contains general information about medical conditions and treatments. The information is not advice and should not be treated as such. This information is not intended to be patient education, does not create any patient-physician relationship, and should not be used as a substitute for professional diagnosis and treatment.

The medical information in this lecture is provided "as is" without any representations or warranties, express or implied. Precision Analytical makes no representations or warranties in relation to the medical information in this presentation.



# Case 5: Anna

Endometriosis and Dysmenorrhea



# Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea

## Chief Complaints:

- Dysmenorrhea
- Heavy bleeding
- High stress
- Headaches

## PMHx

- Stress job x 7 months
- Regular 26-day cycles
- Menses increasingly more painful/heavier in the past 4 months
- Headaches day prior to menses lasting through days 1-2

## Medications

- None

## Physical Exam

- 5'5"; 157 lbs.; BMI 26.1
- BP 116/76 mm Hg
- Pulse 84 bpm

## Pertinent Serum Labs:

- Iron deficiency anemia
- Day 3 Labs: FSH 5.2 mIU/ml, E2 58 pg/ng **(H)**
- Vitamin B12: 345 pg/mL **(L)**
- CRP-hs 3.2 mg/dL **(H)**





# Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea






# Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea





# Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea



## Sex Hormones & Metabolites


TEST	RESULT	UNITS	LUTEAL*	POSTMENOPAUSAL
<b>Progesterone Metabolites (Urine)</b>				
b-Pregnanediol	Within luteal range	1659.0	ng/mg	600 - 2000
a-Pregnanediol	Within luteal range	472.0	ng/mg	200 - 740
<b>Estrogens and Metabolites (Urine)</b>				
Estrone (E1)	Below luteal range	9.53	ng/mg	12 - 26
Estradiol (E2)	Within luteal range	2.50	ng/mg	1.8 - 4.5
Estrinol (E3)	Within luteal range	10.2	ng/mg	5 - 18
2-OH-E1	Below luteal range	1.99	ng/mg	5.1 - 13.1
4-OH-E1	Within luteal range	0.38	ng/mg	0 - 1.8
16-OH-E1	Above luteal range	2.70	ng/mg	0.7 - 2.6
2-Methoxy-E1	Below luteal range	0.87	ng/mg	2.5 - 6.5
2-OH-E2	Within luteal range	0.86	ng/mg	0 - 3.1
4-OH-E2	Within luteal range	0.15	ng/mg	0 - 0.52
Total Estrogen	Below range	29.2	ng/mg	35 - 70
<b>Metabolite Ratios (Urine)</b>				
2-OH / 16-OH-E1 Balance	Below range	0.74	ratio	2.69 - 11.83
2-OH / 4-OH-E1 Balance	Below range	5.24	ratio	5.4 - 12.62
2-Methoxy / 2-OH Balance	Within range	0.44	ratio	0.39 - 0.67
<b>Androgens and Metabolites (Urine)</b>				
<b>Range</b>				
DHEA-S	Within range	55.6	ng/mg	20 - 750
Androsterone	Within range	1318.8	ng/mg	200 - 1650
Etiocholanolone	Within range	994.2	ng/mg	200 - 1000
Testosterone	Within range	12.87	ng/mg	2.3 - 14
5a-DHT	Within range	2.0	ng/mg	0 - 6.6
5a-Androstenediol	Within range	16.2	ng/mg	6 - 30
5b-Androstenediol	Within range	50.0	ng/mg	12 - 75
Epi-Testosterone	Within range	12.6	ng/mg	2.3 - 14

\* The Luteal Range represents the expected premenopausal luteal range, collected menstrual cycle days 19-22 of a 28-day cycle. If your patient noted taking oral progesterone, the reference range represents the expected range on 100 - 200 mg of oral micronized progesterone (OMP). The ranges in the table below represent ranges in other times of the cycle your patient may have collected, such as follicular or ovulatory phases.

ADDITIONAL NORMAL RANGES	FOLLICULAR	OVULATORY	ON ORAL PG
b-Pregnanediol	100 - 300	100 - 300	2000 - 9000
a-Pregnanediol	25 - 100	25 - 100	580 - 3000
Estrone (E1)	4.0 - 12.0	22 - 68	N/A
Estradiol (E2)	1.0 - 2.0	4.0 - 12.0	N/A



# Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea




## Sex Hormones & Metabolites

TEST	RESULT	UNITS	LUTEAL*	POSTMENOPAUSAL
<b>Progesterone Metabolites (Urine)</b>				
b-Pregnanediol	Within luteal range	1659.0	ng/mg	600 - 2000
a-Pregnanediol	Within luteal range	472.0	ng/mg	200 - 740
<b>Estrogens and Metabolites (Urine)</b>				
Estrone (E1)	Below luteal range	9.53	ng/mg	12 - 26
Estradiol (E2)	Within luteal range	2.50	ng/mg	1.8 - 4.5
Estrinol (E3)	Within luteal range	10.2	ng/mg	5 - 18
2-OH-E1	Below luteal range	1.99	ng/mg	5.1 - 13.1
4-OH-E1	Within luteal range	0.38	ng/mg	0 - 1.8
16-OH-E1	Above luteal range	2.70	ng/mg	0.7 - 2.6
2-Methoxy-E1	Below luteal range	0.87	ng/mg	2.5 - 6.5
2-OH-E2	Within luteal range	0.86	ng/mg	0 - 3.1
4-OH-E2	Within luteal range	0.15	ng/mg	0 - 0.52
<b>Total Estrogen</b>	<b>Below range</b>	<b>29.2</b>	<b>ng/mg</b>	<b>35 - 70</b>
<b>Metabolite Ratios (Urine)</b>				
2-OH / 16-OH-E1 Balance	Below range	0.74	ratio	2.69 - 11.83
2-OH / 4-OH-E1 Balance	Below range	5.24	ratio	5.4 - 12.62
2-Methoxy / 2-OH Balance	Within range	0.44	ratio	0.39 - 0.67
<b>Androgens and Metabolites (Urine)</b>				
Range				
DHEA-S	Within range	55.6	ng/mg	20 - 750
Androsterone	Within range	1318.8	ng/mg	200 - 1650
Etiocholanolone	Within range	994.2	ng/mg	200 - 1000
Testosterone	Within range	12.87	ng/mg	2.3 - 14
5a-DHT	Within range	2.0	ng/mg	0 - 6.6
5a-Androstenediol	Within range	16.2	ng/mg	6 - 30
5b-Androstenediol	Within range	50.0	ng/mg	12 - 75
Epi-Testosterone	Within range	12.6	ng/mg	2.3 - 14

\* The Luteal Range represents the expected premenopausal luteal range, collected menstrual cycle days 19-22 of a 28-day cycle. If your patient noted taking oral progesterone, the reference range represents the expected range on 100 - 200 mg of oral micronized progesterone (OMP). The ranges in the table below represent ranges in other times of the cycle your patient may have collected, such as follicular or ovulatory phases.

ADDITIONAL NORMAL RANGES	FOLLICULAR	OVULATORY	ON ORAL PG
b-Pregnanediol	100 - 300	100 - 300	2000 - 9000
a-Pregnanediol	25 - 100	25 - 100	580 - 3000
Estrone (E1)	4.0 - 12.0	22 - 68	N/A
Estradiol (E2)	1.0 - 2.0	4.0 - 12.0	N/A

## Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea




### Organic Acid Tests (OATs)

TEST	RESULT	UNITS	NORMAL RANGE
<b>Nutritional Organic Acids (Urine)</b>			
Vitamin B12 Marker - May be deficient if high			
Methylmalonate (MMA)	High end of range	2.4 ug/mg	0 - 2.5
Vitamin B6 Markers - May be deficient if high			
Xanthurenate	Above range	1.34 ug/mg	0.12 - 1.2
Kynurenate	Above range	4.7 ug/mg	0.8 - 4.5
Biotin Marker - May be deficient if high			
b-Hydroxyisovalerate	Within range	8.5 ug/mg	0 - 12.5
Glutathione Marker - May be deficient if high			
Pyroglutamate	Within range	40.0 ug/mg	28 - 58
Gut Marker - Potential gut putrefaction or dysbiosis if high			
Indican	Above range	103.4 ug/mg	0 - 100
<b>Neuro-Related Markers (Urine)</b>			
Dopamine Metabolite			
Homovanillate (HVA)	Within range	7.8 ug/mg	3 - 11
Norepinephrine/Epinephrine Metabolite			
Vanilmandelate (VMA)	Within range	4.2 ug/mg	2.2 - 5.5
Neuroinflammation Marker			
Quinolinatate	Within range	4.1 ug/mg	0 - 9.6
<b>Additional Markers (Urine)</b>			
Melatonin - Waking			
6-OH-Melatonin-Sulfate	Within range	57.8 ng/mg	10 - 85
Oxidative Stress / DNA Damage			
8-Hydroxy-2'-deoxyguanosine (8-OHdG)	Above range	5.9 ng/mg	0 - 5.2

- Both the xanthurenate and kynurenate are above the range. This may indicate vitamin B6 deficiency. B6 is important for phase 2 methylation (estrogen detox), neurotransmitter synthesis, and other key processes. Tryptophan taken within 72 hours before testing can also raise these markers without indicating a true B6 deficiency.
- The indican is above the range. This can indicate gut dysbiosis. Gut dysbiosis can affect estrogen metabolism, inflammation, and malabsorption of nutrients. Further GI testing may be indicated.
- The 8-hydroxy-2'-deoxyguanosine (8-OHdG) is above the normal range. This can result from oxidative damage to DNA. This may indicate a need for reducing oxidative damage, reactive oxygen species, and improving antioxidant status.

## Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea



### Organic Acid Tests (OATs)

TEST	RESULT	UNITS	NORMAL RANGE
<b>Nutritional Organic Acids (Urine)</b>			
Vitamin B12 Marker - May be deficient if high			
Methylmalonate (MMA)	High end of range	2.4 ug/mg	0 - 2.5
Vitamin B6 Markers - May be deficient if high			
Xanthurenate	Above range	1.34 ug/mg	0.12 - 1.2
Kynurenate	Above range	4.7 ug/mg	0.8 - 4.5
Biotin Marker - May be deficient if high			
b-Hydroxyisovalerate	Within range	8.5 ug/mg	0 - 12.5
Glutathione Marker - May be deficient if high			
Pyroglutamate	Within range	40.0 ue/mg	28 - 58
<b>Gut Marker - Potential gut putrefaction or dysbiosis if high</b>			
Indican	Above range	103.4 ug/mg	0 - 100
<b>Neuro-Related Markers (Urine)</b>			
Dopamine Metabolite			
Homovanillate (HVA)	Within range	7.8 ug/mg	3 - 11
Norepinephrine/Epinephrine Metabolite			
Vanilmandelate (VMA)	Within range	4.2 ug/mg	2.2 - 5.5
Neuroinflammation Marker			
Quinolinatate	Within range	4.1 ug/mg	0 - 9.6
<b>Additional Markers (Urine)</b>			
Melatonin - Waking			
6-OH-Melatonin-Sulfate	Within range	57.8 ng/mg	10 - 85
<b>Oxidative Stress / DNA Damage</b>			
8-Hydroxy-2-deoxyguanosine (8-OHdG)	Above range	5.9 ng/mg	0 - 5.2

- Both the xanthurenate and kynurenate are above the range. This may indicate vitamin B6 deficiency. B6 is important for phase 2 methylation (estrogen detox), neurotransmitter synthesis, and other key processes. Tryptophan taken within 72 hours before testing can also raise these markers without indicating a true B6 deficiency.
- The indican is above the range. This can indicate gut dysbiosis. Gut dysbiosis can affect estrogen metabolism, inflammation, and malabsorption of nutrients. Further GI testing may be indicated.
- The 8-hydroxy-2'-deoxyguanosine (8-OHdG) is above the normal range. This can result from oxidative damage to DNA. This may indicate a need for reducing oxidative damage, reactive oxygen species, and improving antioxidant status.



# Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea

## Goals of Treatment

### DUTCH Test Goals

- Improve 2/16 and 2/4 balance
- Support for low B6
- Address dysbiosis
- Reduce oxidative stress


### Blood Lab Goals

- Improve anemia
- Improve underlying inflammation/elevated CRP

### Lifestyle Goals

- Improve diet and hydration
- Start exercise program

## Treatment Guide Key:



○ HPO Axis Support

○ HPA Axis Support

○ Other Hormone support

○ OATs Support

○ Symptom Support

○ Detox Support

○ Lifestyle Support

○ Other Support

# A Case 5: Anna, a 20-yo Female with Endometriosis & Dysmenorrhea

## Sample Treatment Options:

### Supplements

- Consider **Beef liver** capsules, 2 capsules BID
  - **To address iron-deficiency anemia**
- Consider **herbal pain and inflammation relief supplement** containing rosemary, turmeric, ginger, holy basil, green tea, hu zhang, goldthread, barberry, oregano, skullcap
  - **To support phase 1 estrogen detox, decrease inflammation & oxidative stress, pain relief**
- Consider **Broccoli seed extract blend** that contains glucoraphanin and myrosinase – one capsule daily
  - **To support phase 1 estrogen detox**

### Diet

- Consider dietary intake of healthy fats/proteins daily, increase cruciferous veggie and broccoli sprout intake, add ground flaxseeds daily
  - **To help weight loss efforts, decrease insulin resistance, improve sleep, improve estrogen detoxification**
  - **To support balances E1 to E2 conversion**

### Lifestyle

- Consider 4-5 days a week of moderate intensity exercise – incorporating weight/resistance training
  - **To improve weight loss and help lower estrogens, improve insulin resistance, decrease inflammation and stress**
- Consider meditation and mindfulness program of choice **to optimize cortisol levels, if indicated.**
- Referral Functional GI Health eval
  - **To address for underlying gut dysbiosis and inflammation, which can in turn improve estrogen metabolism and high estrogen levels**

■ DUTCH Dozen

■ Advanced Insights



## Polling Question:

What was the key finding that influenced your treatment plan the most?

- A. Poor phase 1 estrogen detox
- B. Dysbiosis
- C. Oxidative stress
- D. B6 and B12 deficiencies

## **Reflection:**

**How does her slow estrogen clearance in the context of her dysbiosis influence your original treatment plan?**



- 1. How do you know your treatment plan is working?**
- 2. When would you retest?**
- 3. Which DUTCH panel would you order?**

# Thank You!

DUTCH Fest 2026

